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STORAGE AND STABILITY

Store the test kit at 4-30°C with a valid period of 24 months.
Use the test card within 1 hour once the foil pouch is opened.

PRINCIPLE

The test uses anti-SARS-CoV-2 nucleocapsid protein (N protein) monoclonal antibody I conjugated with colloidal gold coated on the sample pad, and another anti-SARS-CoV-2 N protein monoclonal antibody II coated on test line. After the samples have been applied to the test strip, the colloidal gold-labelled anti-SARS-CoV-2 N protein monoclonal antibody I bind with SARS-CoV-2 antigens in sample and form marked antigen-antibody complexes. These complexes move to the test card detection zone by capillary action. Then marked antigen-antibody complexes will be captured on test line by anti-SARS-CoV-2 N protein monoclonal antibody II. The color intensity of each test line increases in proportion to the amount of SARS-CoV-2 antigen in sample.

PRECAUTIONS

1. Always keep out of the reach of children. Small parts of the kit can be a choking hazard.
2. Sample extraction solution is a phosphate buffer contained low concentration of sodium chloride, Tween, hexadecyl trimethyl ammonium bromide and sodium azide. If extraction solution splashes your body or into eyes, please wash with water.

LIMITATIONS

1. False-negative result may occur if the level of antigen in sample is below the detection limit of the test or the sample was collected incorrectly.
2. The test kit cannot differentiate between SARS-CoV and SARS-CoV-2.
3. Clinical diagnosis and treatment cannot be made without consulting with the physician.
4. A negative result, from an individual have symptoms similar to COVID-19 beyond seven days should be treated as negative possibly. If necessary, confirmed with the molecular assay.

PERFORMANCE CHARACTERISTICS

1 Limit of Detection (LoD)

The LoD for saliva was established using heat-inactivated SARS-CoV-2 isolate strain. The strain was spiked with negative human saliva into a series of concentrations. Each level was tested for 20 replicates. The confirmed LoD for saliva was 200 TCID₅₀/mL.

2 Clinical Agreement Study

The clinical performance of One Step Test for SARS-CoV-2 Antigen (Colloidal Gold) (Saliva) was evaluated by testing a total of 216 saliva samples. It is compared to the results of RT-PCR assays. Overall study results were shown in the tables below.

Total		RT-PCR		
		positive	negative	Subtotal
Getein's kit	positive	105	0	105
	negative	1	110	111
	subtotal	106	110	216

Positive percent agreement (Diagnostic sensitivity) = 105/ (105+1) ×100% =99.06% (95% CI: 94.85%-99.83%)

Negative percent agreement (Diagnostic specificity) = 110/ (110+0) × 100%=100.00% (95% CI: 96.63%-100.00%)

Total percent agreement = (105+110)/216× 100%=99.54% (95% CI: 97.42%-99.92%)

3 Analytical Specificity

3.1 Cross-Reactivity & Microbial Interference

Each organism and virus were tested in triplicate in the absence and presence of SARS-CoV-2 respectively. According to the test results, there was no cross-reactivity with the following viruses or organisms.

Viruses or organisms	Concentration
Human coronavirus 229E	1 x 10 ⁵ PFU/mL
Human coronavirus OC43	1 x 10 ⁵ PFU/mL
Human coronavirus NL63	9.87 x 10 ³ PFU/mL
MERS coronavirus	7930 PFU/mL
Adenovirus (e.g. C1 Ad. 71)	1 x 10 ⁵ PFU/mL
Human Metapneumovirus (hMPV)	1 x 10 ⁵ PFU/mL
Parainfluenza virus Type 1	1 x 10 ⁵ PFU/mL
Parainfluenza virus Type 2	1 x 10 ⁵ PFU/mL
Parainfluenza virus Type 3	1 x 10 ⁵ PFU/mL
Parainfluenza virus Type 4a	1 x 10 ⁵ PFU/mL
Influenza A	1 x 10 ⁵ PFU/mL
Influenza B	2.92 x 10 ⁴ PFU/mL
Enterovirus	1 x 10 ⁵ PFU/mL
Respiratory syncytial virus	1 x 10 ⁵ PFU/mL
Rhinovirus	4.17 x 10 ⁵ PFU/mL
Haemophilus influenzae	1 x 10 ⁶ CFU/mL
Streptococcus pneumoniae	1 x 10 ⁶ CFU/mL
Streptococcus pyogenes	1 x 10 ⁶ CFU/mL
Candida albicans	1 x 10 ⁶ CFU/mL
Bordetella pertussis	1 x 10 ⁶ CFU/mL
Mycoplasma pneumoniae	1 x 10 ⁶ CFU/mL
Chlamydia pneumoniae	1 x 10 ⁶ CFU/mL
Legionella pneumophila	1 x 10 ⁶ CFU/mL
Mycobacterium tuberculosis	1 x 10 ⁶ CFU/mL
Pneumocystis jirovecii	1 x 10 ⁶ CFU/mL
Pseudomonas Aeruginosa	1 x 10 ⁶ CFU/mL
Staphylococcus Epidermidis	1 x 10 ⁶ CFU/mL
Streptococcus Salivarius	1 x 10 ⁶ CFU/mL

3.2 Interferences

The potentially interfering substances that may be found in the upper respiratory tract in subjects (including over the counter medications). No false positive or false negative results were seen at the following concentrations.

Potentially Interfering Substances	Concentration
Blood (human)	5%
Mucin	5 mg/mL
Naso GEL (NeilMed)	5% v/v
CVS Nasal Drops (phenylephrine)	15% v/v
Afrin (Oxymetazoline)	15% v/v
CVS Nasal Spray (Cromolyn)	15% v/v
Zicam Cold Remedy	5% v/v
Homeopathic (Alkalol)	10 % v/v
Sore Throat Phenol Spray	15% v/v
Antibacterial, systemic (Tobramycin)	3.3 mg/dL
Antibiotic, nasal ointment (Mupirocin)	0.15 mg/dL
Nasal corticosteroids (Fluticasone)	5% v/v
Anti-viral drugs (Tamiflu-Osetamivir phosphate)	500 mg/dL
Biotin	0.35 mg/dL
Methanol	0.15% w/v
Diphenhydramine	0.0774 mg/dL
Dextromethorphan	0.00156 mg/dL
Dexamethasone	1.2 mg/dL

4.Precision

For repeatability study, the agreement percent of both negative samples and positive samples are 100%. For reproducibility study, the agreement percent of both negative samples and positive samples are 100%.

DESCRIPTION OF SYMBOLS USED

Key to symbols used			
	Manufacturer		Use-by date
	Do not re-use		Date of manufacture
	Consult instructions for use		Batch code
	Temperature limit		In vitro diagnostic medical device
	Contains sufficient for <n> tests		Authorized representative in the European Community
	Keep away from sunlight		Do not use if package is damaged
	Catalogue number		Keep away from rain
	For self-testing		

Getein Biotech, Inc.
Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China
Tel: +86-25-68568508
Fax: +86-25-68568500
E-mail: tech@getein.com.cn overseas@getein.com.cn
Website: en.bio-gp.com.cn

CMC Medical Devices & Drugs S.L.
Add: C/ Horacio Lengo N° 18, CP 29006, Málaga, Spain
Tel: +34951214054

Version: WCG93-DXSX-S-01

Specification (N)	REF
1 T/kit	CG20616
2 T/kit	CG206162
3 T/kit	CG206163
5 T/kit	CG206165
6 T/kit	CG206166
7 T/kit	CG206167
8 T/kit	CG206168
9 T/kit	CG206169
10 T/kit	CG2061610
12 T/kit	CG2061612
15 T/kit	CG2061615
20 T/kit	CG2061620
25 T/kit	CG2061625

Scan this QR Code for Audiovisual Demonstration.

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